

ANALIZA BIAŁEK



Scan for verify



100017890

Laboratory of Analysis

Client: Marcin Sowinski
Order number: 100017890
Name: NIEOZNACZONA
Number of samples: 1 fiołka

Batch number: -----
Receive date: 2026-09-01
Start date: 2026-09-07
End date: 2026-09-09

Sample identification

Substance	Method of taking a sample	Form
Tirzepatide	Próbka dostarczona przez klienta	Biały proszek w szklanej fiołce


Results

Substance	Sample name	Content	Purity
Tirzepatide	NIEOZNACZONA	10,93mg	99%
Endotoxin Level	NIEOZNACZONA	5,83EU/vial	-----

Research comment

The test sample was applied in the amount of 3 μ l to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H₂O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its content, when compared to the reference chromatogram, is 10,93mg. Height the measurement error is $\pm 0,05$ mg. The actual state of the sample sent before the test is presented on the photo below.

Photo of the test sample




Password

LZ34VOM8

Analiza Bialek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zajęc



ANALIZA BIAŁEK




Scan for verify



100017890

Chromatogram



Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	95%	5%
2min	95%	5%
12min	5%	95%
15min	5%	95%

Password

LZ34VOM8

Analiza Białek sp. z o.o.

KRS: 0000966168

NIP: 8992922970

REGON: 521751875

Headquarters and Laboratory:

ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:

Analysis and research specialist
mgr Arkadiusz Zając